

## Atlantis® suprastructures compatibility chart - extended assortment



Manufacturer/Company	Implant system	Screw Order REF number from Manufacturer/Company*
Alphatech	<b>IMPLANT LEVEL</b>	
	D3.4, D3.8, D4.3 Implant Level	188705 (5 screws)
	D5, D6, D7 Implant Level	188705 (5 screws)
	<b>ABUTMENT LEVEL</b>	
	Abutment Level	188704 (5 screws)
Alpha Bio	<b>IMPLANT LEVEL</b>	
	Internal Hex	STLAS
	<b>ABUTMENT LEVEL</b>	
	TSA Abutment	SF
Anthoghyr	<b>ABUTMENT LEVEL</b>	
	4.0 Multi-Unit Abutment Narrow	MU140
	4.8 Multi-Unit Abutment	MU140
	Conical Abutment D4.0 (Old)	OPAC141
	Conical Abutment D4.8 (Old)	OPAC141
Bego	<b>IMPLANT LEVEL</b>	
	Semados SC/SCX/RS/R SX/RI 3.25, 3.75	57031
	Semados SC/SCX/RS/R SX/RI 4.1	57031
	Semados SC/SCX/RS/R SX/RI 4.5	57031
	Semados SC/SCX/RS/R SX/RI 5.5	57031
	Semados S/RI D3.25/D3.75	57031
	Semados S/RI D4.1, D4.5, D5.5	57031
	Semados Mini	56983
	<b>ABUTMENT LEVEL</b>	
	Semados Sub-Tec MultiPlus	57534
BioComp	<b>IMPLANT LEVEL</b>	
	D4 Implant Level	BC-1819-225220
	D5.2 Implant Level	BC-1819-225220
	<b>ABUTMENT LEVEL</b>	
	D4/D5.2 Abutment Level	BC-1824-282015
Biotech	<b>ABUTMENT LEVEL</b>	
	Bis D3/D3.4/D3.6/D3.9/D4.4 Abutment Level	E007
	Conical Abutment D4.9	EAVGC
BT	<b>IMPLANT LEVEL</b>	
	Lock D4.5	370 BTIVPT
BTI	<b>IMPLANT LEVEL</b>	
	Tiny D3.5	TTTH
	Interna Universal D4.1	INTTUH
	Interna Wide D5.5	INTTUH
	<b>ABUTMENT LEVEL</b>	
	Multi-IM Abutment	TTMI
Bredent	<b>IMPLANT LEVEL</b>	
	SKY D4.1 Implant Level	SKY PS22 (6X)
	<b>ABUTMENT LEVEL</b>	
	SKY Fast/Fixed Abutment Level	SKYFFSPK (6X)
	SKY uni.cone abutment	SKYFFSPK (6X)

\* Note: We recommend to check Screw REF number with the manufacturer.

All trademarks, company names and implant designs are the property of their respective owners.

Please turn over →

Manufacturer/Company	Implant system	Screw Order REF number from Manufacturer/Company*
Dentaurum	<b>IMPLANT LEVEL</b>	
	TioLogic S Implant Level	386-910-00
	TioLogic M Implant Level	386-910-00
	TioLogic L Implant Level	386-910-00
	TioLox D3.5, D4.5 Implant Level	330-790-00 (5mm)
	<b>ABUTMENT LEVEL</b>	
	TioLogic Abutment Level	386-912-00
Dentegris	<b>IMPLANT LEVEL</b>	
	D3.3 Implant Level	ZS 01
	D3.75 Implant Level	ZS 01
	D4.5/D4.75 Implant Level	ZS 01
	D5.5 Implant Level	ZS 03
DIO	<b>IMPLANT LEVEL</b>	
	SM-Submerged D3.8/4.1 Implant Level Narrow	SSC 1808
	SM-Submerged D4.5/D5.0/D5.3 Implant Level Regular/Wide	SSC 2008
	SM-Internal D4.8 Implant Level Regular/Wide	STN 2008
Implant Direct	<b>IMPLANT LEVEL</b>	
	Legacy D3	1000-04
	Legacy D3.5 / 4.5 / 5.7	1000-01
	S-plant D3.7, D4.7, D5.7	1000-01
	S-InDirect D5.0	1000-01
	<b>ABUTMENT LEVEL</b>	
	Overdenture Screw-Retained Abutment	1000-05
Imtec	<b>IMPLANT LEVEL</b>	
	Implant-Level D3.5/D4.3	5052HX
	Implant-Level D5.1	5052HX
Medical Instinct	<b>IMPLANT LEVEL</b>	
	BoneTrust D3, D3.4 Implant Level	160-100004
	BoneTrust D4, D5 Implant Level	160-100004
	<b>ABUTMENT LEVEL</b>	
	D3.4 /D4 Direct Abutment	165-344055
Megagen	<b>IMPLANT LEVEL</b>	
	Exfeel Internal Implant Level	EXIMS100
	Rescue Internal Implant Level	ARAS250L
	Rescue External D5 Implant Level	LULA300
	EZ Plus Internal Implant Level Small	EPAS100
	EZ Plus Internal Implant Level Regular/Wide	AMAS100
	<b>ABUTMENT LEVEL</b>	
	Exfeel Internal Abutment Level	IRCS200
Medentis	<b>IMPLANT LEVEL</b>	
	ICX Templant	C-011-000001-01
Neobiotech	<b>IMPLANT LEVEL</b>	
	IS-II active/IS-II/IS Narrow	ISCS20
	IS-II active/IS-II/IS Regular/Wide	ISCS20
	<b>ABUTMENT LEVEL</b>	
	IS Multi Unit Abutment	MUAS 14
Neoss	<b>IMPLANT LEVEL</b>	
	Narrow D3.25 Implant Level	31240
	D3.5/D4.0/D4.5/D5.0/D5.5 Implant Level	31150
	<b>ABUTMENT LEVEL</b>	
	D4.1 Access Abutment	31229

\* Note: We recommend to check Screw REF number with the manufacturer.

Please turn over 

All trademarks, company names and implant designs are the property of their respective owners.



Manufacturer/Company	Implant system	Screw Order REF number from Manufacturer/Company*
Pitt Easy	<b>IMPLANT LEVEL</b>	
	D3.25, D3.75, D4, D4.9, D6.5	90276
Platon	<b>IMPLANT LEVEL</b>	
	PF D4.4 Implant Level	QPB-0001
	PF D4.8 Implant Level	QPB-W001
	Pro PF D4.1, D5.0 Implant Level	QPB-0001
	Eight-Lobe NP AM	QRB-0044
	Eight-Lobe RP AM	QRB-0048
	<b>ABUTMENT LEVEL</b>	
	PF D4.1/D4.4 Abutment Level	QPA-0005/QPB-0001
	PF D4.8 Abutment Level	QPA-W005/QPB-W001
POI	<b>IMPLANT LEVEL</b>	
	EX D3.7 Implant Level	CD012600
	EX D4.2 Implant Level	CD012650
	EX D4.7 Implant Level	CD012700
	EX D5.2 Implant Level	CD012750
Prowital	<b>IMPLANT LEVEL</b>	
	Wieland D3.5, D4.3, D5 Implant Level	12.131.600
Renew Biocare	<b>ABUTMENT LEVEL</b>	
	UMA Abutment	11-000052
SIC	<b>IMPLANT LEVEL</b>	
	D3.3, D4.2, D5.3	936658
	<b>ABUTMENT LEVEL</b>	
	Multi-Unit Abutment	936540
	Mini Multi-Unit Abutment	936295
Simpler	<b>ABUTMENT LEVEL</b>	
	Bar Abutment	AS 21001
Serf	<b>IMPLANT LEVEL</b>	
	EVL Plus D3.3	VPEG 3,3-DP30700050
	EVL Plus D4, D5	VPEG 4 ET 5-DP30700070
Schuetz	<b>IMPLANT LEVEL</b>	
	Dental Impla D3.3, D4.2, D5.3	636649
Sweden Martina	<b>IMPLANT LEVEL</b>	
	Global D3.8, D4.3, D4.8, D5.5	U-VMOU1.8
	Premium/Kohno D3.3, D3.8	VM2-180
	Premium/Kohno D4.25, D5	VM2-200
	<b>ABUTMENT LEVEL</b>	
	PAD Abutment Level	PAD-VP-140
Tekka	<b>IMPLANT LEVEL</b>	
	Progress D4.5, D5, D6 Implant Level	DVPQP
	Krestal D3.6 Implant Level	DVPQP
	Monobloc D3.4, D4, D4.5, D5	DVPIP
	Quattro	DVPAOFP
Thommen Medical	<b>IMPLANT LEVEL</b>	
	SPI D3.5 Implant Level	4.03.500
	SPI D4, D4.5, D5, D6 Implant Level	4.03.501
	<b>ABUTMENT LEVEL</b>	
	SPI Abutment Level	4.03.506
	SPI VARIO multi Abutment	4.03.507

\*Note: We recommend to check Screw REF number with the manufacturer.

Note: Atlantis suprastructures assortment can vary from country to country. Please check Atlantis WebOrder for availability.  
All trademarks, company names and implant designs are the property of their respective owners.